


enigma



CONVERSION CHART



DILLON Co., INC.

LEACH & DILLON PRODUCTS


Canyon State
Dental Supply
1-800-354-3684
1-800-661-2044


schottlander
committed to quality

K27	49.2	10.9	9.1	T99	48.0	11.1	9.0	45H	49.5	11.25	8.85	A27	48.5	10.9	8.9
K27	49.2	10.9	9.1	T98	50.5	10.8	9.3								
K41	42.7	9.8	8.2	T36	43.7	10.1	8.1	42C	42.00	9.5	7.5	A41	41.8	9.9	8.1
K42	44.3	9.6	8	T66	44.6	10.2	8.3	43F	45.5	9.6	8.8	A42	44.0	9.9	7.8
K42	44.3	9.6	8	T53	44.6	9.2	8.1								
K42	44.3	9.6	8	T67	45.1	10.0	8.4								
S11	41.5	8.6	8.1	Z51	44.6	8.8	8.3	13E	44.0	9.3	8.1	A11	41.3	8.4	8.0
S12	45.1	8.5	8.6	Z61	45.5	8.7	8.2	13E	44.0	9.3	8.1	A12	45.0	8.8	8.4
S12	45.1	8.5	8.6									A64	46.9	8.6	8.6
S13	45.1	8.5	8.6					13E	44.0	9.3	8.1				
S13	43.1	9.4	8.4	Z51	44.6	8.8	8.3	13D	44.0	9.5	7.75	A13	43.5	9.8	8.3
S14	46.6	10.7	8.9	Z74	46.7	9.6	8.8	12E	45.5	10.6	8.5	A14	46.3	10.8	8.8
S14	46.6	10.7	8.9	Z84	46.7	9.4	8.7								
S15	50.3	10.7	9.4	N/A	N/A	N/A	N/A	12G	49	10.6	9.15	A15	49.1	10.8	9.0
S15	50.3	10.7	9.4	N/A	N/A	N/A	N/A	11G	48	11.1	8.75	A68	49.4	9.7	9.2
S17	52	12.8	9.8	Z97	52.2	10.2	9.5	11H	50.50	11.60	9.35	A17	52.0	13.0	9.8
S66	46.4	10.3	8.8	Z85	47.5	9.7	8.6	12F	47.0	10.25	8.5	A66	46.0	10.3	8.7
O31	43	9.3	7.9	O13	41.7	9.5	7.9	62D	43.0	10.1	8.1	A31	42.7	9.4	7.8
O31	43	9.3	7.9	O25	42.8	9.9	8.1								
O31	43	9.3	7.9	O43	43.6	9.4	8.2								
O32	45.1	9.5	8.4	O34	43.3	9.7	8.0	62E	44.0	10.0	8.0	A32	44.7	9.5	8.3
O32	45.1	9.5	8.4	O35	43.2	9.8	8.0								
O34	48.4	9.8	8.8	O86	47.1	9.9	8.9	62G	48.5	11.9	9.2	A34	48.0	9.9	8.6
O36	46.9	10.3	8.7	O86	47.1	9.9	8.9	65G	47.5	11.5	9.0	A36	46.4	10.3	8.6
O36	46.9	10.3	8.7	O98	48.1	11.8	9.0								
O36	46.9	10.3	8.7	O97	49.2	11.5	9.0								
O54	47.1	9.5	8.9	O86	47.1	9.9	8.9	65G	47.5	11.5	9.0	A54	46.3	9.6	8.6
N/A	N/A	N/A	N/A	O99	52.5	12.5	9.8								

* Length of anteriors are measured flat

COMBINATION TABLES

ANTERIORS		POSTERIORS
Upper	Lower	
S11	L3	P3
S12	L5	P3
S13	L5	P3/P5
S14	L6	P5
S15	L8	P4
S17	L9	P6
S66	L5	P5
K21	L3	P3
K22	L3	P3/P5
K24B	L4	P2
K25	L8	P4
K27	L8	P4
K41	L3	P3/P5
K42	L5	P3/P5
O31	L3	P3
O32	L3	P3/P5
O34	L5	P5
O36	L7	P5
O54	L5	P5
D36	L6	P3/P5
D56	L5	P5
D66	L5	P5
D67	L5	P5
D76	L7	P5
D77	L7	P4/P5
D88	L7	P4/P5
D99	L8	P4

SHADE CONVERSION

Vita®	Bioform® IPN	Bioblend® IPN	Ivoclar®	Chromascope
A1	B 91		1A/120	
A2	B 53	104	2A/130	
A3	B 54	108	1C/140	
A3.5	B 83	116	1E/230	
A4	B 84	118	3E/340	
B1	B 59	100	1A/120	
B2	B 52		1A/120	
B3	B 55		2B/210	
B4	B 56	110	5B/320	
C1	B 62	102	1A/120	
C2	B 94		4B/430	
C3	B 95		6D/510	
C4	B 96		4C/520	
D2	B 92		2A/130	
D3	B 93		4A/410	
D4	B 69	113	6C/440	

enigma is also available in a Hollywood Shade - ask for Shade HB0

DILLON CO., INC.



LEACH & DILLON PRODUCTS

161 COMSTOCK PKWY
CRANSTON, RI 02921
401.464.5850 • 800.535.2633 • 401.464.5862 FX
WWW.LEACHDILLON.COM